

## TREPOR SURGICAL TAPE

*Hypoallergenic*

Non-woven adhesive plaster, breathable, radio transparent, hypoallergenic. Suitable for the fixation of medical devices.



Size	Packaging
1.25 cm X 9.1 m	24 pcs/box
2.50 cm X 9.1 m	12 pcs/box
5.00 cm X 9.1 m	6 pcs/box

## TREPOR SILK TAPE

*Hypoallergenic*

Acetate adhesive silky appearance surgical plaster, breathable, strong and non-elastic, hypoallergenic, easy to tear in both directions by hands, particularly suitable on sensitive skins.



Size	Packaging
2.50 cm X 9.1 m	12 pcs/box
5.00 cm X 9.1 m	6 pcs/box

## TREPOR ZINC OXIDE TAPE

Acetate adhesive silky appearance surgical plaster, breathable, strong and non-elastic, hypoallergenic, easy to tear in both directions by hands, particularly suitable on sensitive skins.



Size	Packaging
2.50 cm X 4.6 m	12 pcs/box
5.00 cm X 4.6 m	6 pcs/box

## TREEFIX+ NON-WOVEN SURGICAL ROLL

*Dressing Retention Tape*

It is very soft, skin respectful and highly hypoallergenic. Allows a comfortable fit of the plaster to the skin and to the other particularly irregular anatomical parts of the body, ensuring the highest moisture permeability. The adhesive is not weakened by water nor secretions. Precut blue squared printed paper protects the adhesive side. Designed for an easy fixation of primary dressings.



Size	Packaging
2.5 cm X 10 m	1 pc/box   320 pcs/carton
5.0 cm X 10 m	1 pc/box   160 pcs/carton
10.0 cm X 10 m	1 pc/box   80 pcs/carton

## FLAT DRESSING

*Polyurethane Adhesive Film*

For a secondary dressing for wounds, widely used in medical preparation fixation, cannula fixation, plaster fixation, fixation of tubing/catheters, etc.

- Ultra-thin (0.01 mm)
- Virtually invisible
- Breathable & Comfortable
- Waterproof


Size	Packaging
10 cm X 10 m	1 roll/box
10 cm X 30 cm (Sterile)	1 pc/pack   100 packs/box



## PARAFFIN GAUZE DRESSING

Paraffin Gauze is soothing and low-adherent and allows the wound to drain freely into an absorbent secondary dressing.

Size	Packaging
10 cm X 10 cm	1 pc/pack   10 pcs/box



## SELF-ADHESIVE BANDAGE

- A soft, stretchy and flexible tape
- Sticks to itself, but will not adhere to skin, hair, clothing or another surface
- Comfortable to wear
- Easy to apply



Size	Packaging
2.5 cm x 4.5 m	1 roll/pack
5.0 cm x 4.5 m	1 roll/pack
7.5 cm x 4.5 m	1 roll/pack
10.0 cm x 4.5 m	1 roll/pack

## CREPE BANDAGE – COTTON

- Suitable for retention bandaging, securing dressings and minor muscle injuries
- Material: 100% cotton
- Colour: Natural white

Size	Packaging
5.0 cm x 4.5 m	12 rolls/pack
7.5 cm x 4.5 m	12 rolls/pack
10.0 cm x 4.5 m	12 rolls/pack
15.0 cm x 4.5 m	12 rolls/pack



## CREPE BANDAGE – HIGH ELASTIC

- Material: Cotton, spandex
- Colour: Brown

Size	Packaging
5.0 cm x 4.5 m	12 rolls/pack
7.5 cm x 4.5 m	12 rolls/pack
10.0 cm x 4.5 m	12 rolls/pack
15.0 cm x 4.5 m	12 rolls/pack



## NITRILE EXAMINATION GLOVE, Powder-Free

<b>Design</b>	Ambidextrous, Beaded Cuff	<b>Protein Level</b>	N/A	
<b>Surface</b>	Palm Textured/ Finger Textured	<b>Powder Content</b>	Max 2.0 mg/glove	
<b>Material</b>	Nitrile Butadiene Rubber	<b>Colour</b>	Blue	
<b>Force at Break</b>	Before Ageing: 9.0-11.0 N After Ageing: 6.0-10.0 N	<b>Elongation</b>	Before Ageing: 650-750 % After Ageing: 500-700 %	
<b>Tensile Strength</b>	Before Ageing: 18-22 MPa After Ageing: 14-18 MPa	<b>Performance Standard</b>	ASTM D-6319, EN 455-Part 1,2&3	
<b>Donning Aid</b>	Chlorinated/ Polymer	<b>Shelf-life</b>	3 Years	
<b>Sizes</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>
<b>Width (mm)</b>	83 ± 6	95 ± 6	106 ± 6	115 ± 6
<b>Length (mm)</b>	Min. 240			
<b>Weight (g)</b>	3.0 ± 0.5	3.5 ± 0.5	4.0 ± 0.5	4.5 ± 0.5
<b>Packaging</b>	100 pcs/box   10 boxes/carton			



## LATEX EXAMINATION GLOVE, Powder-Free

<b>Design</b>	Ambidextrous, Beaded Cuff	<b>Protein Level</b>	Max. 50 Micrograms/ dm2		
<b>Surface</b>	Palm Textured/ Finger Textured/ Smooth	<b>Powder Content</b>	Max 2.0 mg/glove		
<b>Material</b>	Natural Rubber Latex	<b>Colour</b>	Blue		
<b>Force at Break</b>	Before Ageing: 6.0 - 9.0 N After Ageing: 6.0 - 9.0 N	<b>Elongation</b>	Before Ageing: 650 - 800 % After Ageing: 500 - 750 %		
<b>Tensile Strength</b>	Before Ageing: 18 - 22 MPa After Ageing: 14 - 18 MPa	<b>Performance Standard</b>	ASTM D-6319, EN 455-Part 1,2&3		
<b>Donning Aid</b>	Chlorinated/ Polymer	<b>Shelf-life</b>	3 Years		
<b>Sizes</b>	<b>XS</b>	<b>S</b>	<b>M</b>	<b>L</b>	<b>XL</b>
<b>Width (mm)</b>	77 ± 3	83 ± 3	95 ± 3	106 ± 3	112 ± 3
<b>Length (mm)</b>	Min. 240				
<b>Weight (g)</b>	4.5 ± 0.5	5.0 ± 0.5	5.5 ± 0.5	6.0 ± 0.5	6.5 ± 0.5
<b>Packaging</b>	100 pcs/box   10 boxes/carton				



## ALCOHOL PREP PAD

- One pad saturated with 70% Isopropyl Alcohol
- For professional and hospital use
- For disinfection, topical antiseptic use
- For external use only

Type	Size	Packaging
2ply, Non-woven	30x30mm (65x30mm unfolded)	1pc/Aluminium bag <ul style="list-style-type: none"> <li>• 100 pcs/box</li> <li>• 100 boxes/carton</li> </ul>

