

WHEEL CHAIR

- Heavy duty, capacity 130kg
- Extra Large Sizes
- Powder Coated Steel
- Nett Weight 23kg



QUADCANE

WIDE BASE

- Aluminium Framework with Rubber Tip
- Adjustable Height
- Self-Standing



RIPPLE MATTRESS

BUBBLE TYPE

- With Alternating Pressure Pump



IV DRIP STAND



SHARP WASTE BIN

PAIL & COVER

- Litres: 2.5 / 5 / 10 / 20



CLINICAL WASTE BIN

PEDAL OPERATED

- Litres: 10 / 20



HOSPITAL CARE BED

MANUAL 2 CRANKS

- Dimension: L2100 x W900 x H400 mm
- Backrest Tiling: 80° ± 5°
- Knee Board Tilting: 40° ± 5°



HOSPITAL CARE BED

MANUAL 1 CRANK

- Dimension: L2100 x W900 x H400 mm
- Backrest Tiling: 80° ± 5°



OXYGEN CONCENTRATOR

ANGELBISS

- A touch operates multi-functional units
- Large 6" LED Display all functional data
- Timer and accumulated total running hours

Model No	ANGEL-5A	Nebulizer	No
Max O ₂ Flow	5L/min	Humidified Bottle	1
O ₂ Purity	93% ± 3%	Outlet Tube	1
O ₂ Pressure	≤ 0.7 bar	Cannula	1
Noise	≤ 43 dB(A)	Voltages	AC 220-240 V
O ₂ 00.0% Display	Yes	Watts	350W
Timer	Yes	Sizes	380*270*580 mm
Alarm	Yes	Nett Weight	16.5 KG



OXYGEN CYLINDER

3 Liter (0.4 m³)

Aluminium Oxygen Cylinder



5 Liter (0.7 m³)

Aluminium Oxygen Cylinder



10 Liter (1.4 m³)

Aluminium Oxygen Cylinder



10 Liter (1.4 m³)

Steel Oxygen Cylinder

