

# GoChek2 Blood Glucose Monitoring System



### Features

- Meal markers
- Accurate result
- Ketone warning
- Tiny sample size 0.5µl
- Strip ejector
- Auto coding
- High/Low alerts
- Memory 500 readings
- Large display
- **2 units X CR 2032 3.0V** coin cell battery
- Measurement Range : **0.6 - 33.3 mmol/L**
- Haematocrit range : **10% - 70%**

Descriptions	Packaging
Blood Glucose Meter	1 unit
Strip	2 x 25 vials



SCAN ME!

# URIGHT Blood Glucose Monitoring System

TD-4279



### Features

- Ketone warning
- Strip eject function
- Before/after meal (AC/PC) presentable
- USB or Bluetooth data transmission
- Large and backlit display
- Average function 7, 14, 21, 28, 60, 90 days
- **Battery AAA x 2 pcs**
- Measurement Range : **0.55 - 38.89 mmol/L**
- Haematocrit range : **0% - 70%**

Descriptions	Packaging
Blood Glucose Meter	1 unit
Strip	50 pcs



SCAN ME!

# Multi-functional Monitoring System

TD-4289



### Features

- Strip guiding light
- Automatic strip identification system
- LCD with backlight

### 5 types of tests:

- Blood Glucose
- Ketone
- Total Cholesterol
- Uric Acid
- Lactate



SCAN ME!



Blood Glucose



Ketone



Uric Acid




Total Cholesterol



Lactate

# SAFETY LANCET

TD-5052

Gauge	28G	30G
Lancet Colour	Orange 	Purple 
Penetration Depth	1.8mm	1.6mm
Blood volume	5 – 30 µL	3 – 20 µL
Packaging	100 pcs/box 60 boxes/carton	100 pcs/box 60 boxes/carton



SCAN ME!

### Features

- Easy to use
- Single-use only
- Lancet automatically retraction

# BLOOD PRESSURE MONITOR

JPD-HA300



SCAN ME!



Average of last 3 readings

- Accessories
- STORAGE BAG X1
  - AAA BATTERIES X4
  - USER MANUAL X1
  - UNIVERSAL CUFF X1
  - QUICK GUIDE X1

- Heartbeat
- Memory
- Arrhythmia Reminder
- Average
- Voice

# INSULIN PEN NEEDLE

- Compatible with most pens
- Acquired CE, FDA & GMP
- Sterile, Non-Toxic, Non-Pyrogenic

Sizes	Packaging
31G x 5 mm, Purple	100 pcs/box
32G x 4 mm, Green	100 pcs/box
32G x 6 mm, Blue	100 pcs/box




# PORTABLE COMPRESSOR NEBULIZER

TD-7014



SCAN ME!

**DUAL POWER OPTION**



Battery operated



Power adapter

**DC MICROPUMP TECHNOLOGY**



Low-noise



Handheld

**EFFICIENT**

- Fine particle size and effective dose delivery
- Dual filter

Specifications	
Particle size	: Approx. 3 microns
Nebulization rate	: ≥ 0.25 mL/min
Operating pressure range	: 35 - 55 kPa (5 - 8 Psi)
Medication capacity	: 6 ml



# INFECTIOUS DISEASE

## Respiratory Tract Disease



**Fast**  
Test result in 10-20 minutes




**Flexible**  
Able to test different types of patient specimens



**Reliable**  
High sensitivity, specificity and reproductivity



**Intelligent**  
Allows simultaneous screening of different targets



Descriptions	Type	Specimen	Packaging	Flyer
Influenza A & B antigen	Cassette	NS/ NA/ TS	25 test/box	
Influenza A&B + RSV + Adeno	Cassette	NS	25 test/box	

Note: **NS** – Nasal Swab, **NA** – Nasal Aspirate, **TS** – Throat Swab

# INFECTIOUS DISEASE

## Tropical Disease



**Fast**  
Test result in 10-20 minutes




**Flexible**  
Able to test different types of patient specimens



**Reliable**  
High sensitivity, specificity and reproductivity



**Intelligent**  
Allows simultaneous screening of different targets



Descriptions	Type	Specimen	Packaging	Flyer
Dengue NS1 antigen	Cassette	W/S/P	25 test/box	
Dengue IgG/IgM antibodies	Cassette	W/S/P	25 test/box	
Dengue NS1 + IgG/IgM	Cassette	W/S/P	25 test/box	
	Cassette	W/S/P	10 test/box	

Note: **W** – Whole blood, **S** – Serum, **P** – Plasma



## WOMEN'S HEALTH



**Fast**  
Test result in 3 minutes



**Flexible**  
Screening for the presence of hCG in urine in about a week after conception



**Highly Informative**  
No cross-reactivity  
● 300mIU/mL of hLH  
● 300mIU/mL of hFSH  
● 1000mIU/mL of hTSH



Descriptions	Type	Specimen	Packaging	Flyer
hCG Pregnancy 25 mIU/ml	Strip, 4mm	Urine	100 test/box	
hCG Pregnancy 25 mIU/ml	Strip, 5mm	Urine	100 test/box	
hCG Pregnancy 25 mIU/ml	Cassette	Urine	40 test/box	

## INFECTIOUS DISEASE

### Blood Banking



**Fast**  
Test result in 16 minutes



**Accurate**  
High sensitivity, 2NCU/ml



**Reliable**  
Simple and safe to use



**Serology**  
Detection of antibodies to HIV1/2 in human serum, plasma or whole blood



Descriptions	Type	Specimen	Packaging	Flyer
Rapid One Step HIV 1/2 (WHO Pre-Qualified IVD Product List)	Cassette	W/S/P	20 tests/box (Individual pack)	

Note: **W** – Whole blood, **S** – Serum, **P** – Plasma


## DRUG OF ABUSE

Specimen: Urine

Descriptions	Parameter	Packaging	Flyer
5 in 1	AMP, Amphetamine BZO, Benzodiazepine MET, Methamphetamine MOP, Morphine THC, Marijuana	25 test/box	
3 in 1	MET, Methamphetamine MOP, Morphine THC, Marijuana	25 test/box	
Single		50 strip/box	




## ALLTEST URINALYSIS REAGENT STRIP

Descriptions	Parameter	Packaging	Flyer
2 in 1	GLU, Glucose PRO, Protein	100 strips/bottle	
10 in 1	GLU, Glucose PRO, Protein pH, Alkaline LEU, Leukocytes NIT, Nitrite URO, Urobilinogen BLO, Blood SG, Specific Gravity KET, Ketone. BIL. Bilirubin	100 strips/bottle	



## ALLTEST EXCHEK™ ACE URINE ANALYZER

AUR-200

Parameter	Specification	Flyer
Principle	Photoelectric Colorimetry	
Testing speed	60 test/hour	
Data memory	10,000 results	
Display	3.5" LCD touch screen	
Printer	Built-in thermal printer Detachable	
Connectivity	RS232 port, USB ports	

